

Irriguide Lab Times

Irriguide Platinum Full Design Service

	Standard	Rush
Tooth supported: 1-3 implant sites		
First review	1 day	same day
Manufacturing	4 days	2 days
Tooth supported: 4+ implant sites		
First review	2 days	1 day
Manufacturing	4 days	2 days
Bone or mucosa supported: 1-2 implant sites		
First review	2 days	1 day
Manufacturing	5 days	3 days
Bone or mucosa supported 3+ implant sites		
First review	3 days	2 day
Manufacturing	5 days	3 days
Provisionals (add time to guides)		
Full Contour crown	3 days	2 days
Tissue Former	3 days	2 days
Acrylic shell with positioner	3 days	2 days
Essix retainer with pontic	3 days	2 days
Essix retainer without pontic	1 day	1 day
Printed wax up model	1 day	1 day
CBCT radiologist interpretation	7 days	7 days

Irriguide Gold Partial Design Service

Tooth supported		
Manufacturing	4 days	2 days
Tooth supported 4+ implant sites		
Manufacturing	4 days	2 days
Bone or mucosa supported 1-2 implant sites		
Manufacturing	5 days	3 days
Bone or mucosa supported 3+ implant sites		
Manufacturing	5 days	3 days

All required digital files must be uploaded for a case to be considered submitted

Cases submitted before 12:00 Central Time are considered to be submitted that business day

Cases submitted after 12:00 Central Time are considered to be submitted the next business day

